


# Small Headless Screws

CANNULATED DUAL THREAD IMPLANTS




Hand sharpened. Feel the difference.

Inquire about these other precision engineered extremity solutions from Vilex:



1st MPJ Fusion Plate



HammerFUZE™  
Hammertoe Compression  
System



Nitinex™ Memory  
Compression Implant  
and One Cut Blade



Titanex™ Small  
Cannulated Screw  
System

**VILEX**



## Implant sizing

2.0mm, 2.5mm, 3.0mm and 4.2mm implants in lengths to meet most clinical needs

2.0mm	2.5mm	3.0mm	4.2mm
Cannulation = 0.9mm Length = xxT Driver = 1.6mm Hex Drill = Z402-100-17	Cannulation = 1.1mm Length = xxT Driver = 2.0mm Hex Drill = Z402-100-17	Cannulation = 1.1mm Length = xxT Driver = 2.5mm Hex Drill = Z402-100-22	Cannulation = 1.4mm Length = xxT Driver = 2.5mm Hex Drill = Z402-100-27
F20-06T-09	F25-10T-11	F30-10T-11	F42-16T-14
F20-07T-09	F25-12T-11	F30-12T-11	F42-18T-14
F20-08T-09	F25-14T-11	F30-14T-11	F42-20T-14
F20-09T-09	F25-16T-11	F30-16T-11	F42-22T-14
F20-10T-09	F25-18T-11	F30-18T-11	F42-24T-14
F20-11T-09	F25-20T-11	F30-20T-11	F42-26T-14
F20-12T-09	F25-22T-11	F30-22T-11	F42-28T-14
F20-13T-09	F25-24T-11	F30-24T-11	F42-30T-14
F20-14T-09	F25-26T-11	F30-26T-11	F42-32T-14
F20-15T-09	F25-28T-11	F30-28T-11	F42-34T-14
F20-16T-09	F25-30T-11	F30-30T-11	F42-36T-14
F20-17T-09	F25-36T-11	F30-32T-11	F42-38T-14
F20-18T-09		F30-34T-11	F42-40T-14
F20-19T-09		F30-36T-11	F42-45T-14
F20-20T-09		F30-38T-11	F42-50T-14
F20-21T-09		F30-40T-11	F42-55T-14
F20-22T-09		F30-42T-11	F42-60T-14
F20-23T-09		F30-44T-11	
F20-24T-09			
F20-25T-09			
F20-26T-09			
F20-27T-09			
F20-28T-09			
F20-29T-09			
F20-30T-09			
F20-32T-09			
F20-34T-09			
F20-36T-09			
F20-38T-09			
F20-40T-09			
F20-45T-09			
F20-50T-09			



## System Advantages

- Aggressive thread depth enhances pull-out strength and purchase
- Large thread pitch differential maximizes compression
- Conically shaped head and fine threads to engage near cortex
- Hand-sharpened cutting tips, self-drilling, and self-tapping
- Titanium implants offer strength and biocompatibility
- 3.0mm implants available with short or long threads in select sizes
- 2.5mm Stainless Steel implants available in 10mm to 30mm lengths



ADVANCE WITH US >>>

111 Moffitt Street  
McMinnville, TN 37110

p: 800.521.5002  
f: 866.606.4911

info@vilex.com  
www.vilex.com